

Release Notes

Viedoc 4.9

Available for users on
2014-10-10

Release Notes

Viedoc 4.9

New functionality



List of new functionality that has been added to this release.

Designer

- It is now possible to generate a PDF with empty form designs (empty CRFs) with or without annotations. This functionality allows for all forms in the design or an individual form to be exported.

Admin

- It is now possible for study administrators (Study and Site managers) to send reset password emails to users that need to reset their password and have forgotten their security question.
- The listing of study designs in ODM study data import now shows the version and revision number for each available design.
- A user invited to multiple study sites with one or multiple roles will only receive one email message listing all the study sites and roles that the user has been assigned.

Release Notes

Viedoc 4.9

Updated functionality and bug fixes



List of new functionality that has been updated in this release.

- Bugfix: Not possible to enter numbers using the number keys (shift+number) for users with an AZERTY keyboard layout.
- Bugfix: Event status incorrectly shown as 'Ready' for events where all forms have been filled but at least one form contains errors.
- Bugfix: Column headers not enclosed in quotes when column label contains a comma or quotes.
- Bugfix: ViedocMe questionnaires not shown when subject activated with Turkish language

Release Notes

Viedoc 4.9

Known limitations



List of limitations that have been found with this and earlier releases of Viedoc, and will be changed in future releases.

- The export preview feature is heavy to load and might cause the application to feel unresponsive
- Only one (selected) language can be imported from each CDISC ODM file
- Only one (selected) design can be imported from each CDISC ODM file
- When the study workflow is saved, the language in the design settings are lost and have to be reset
- Visibility conditions (entry conditions) cannot be set for activities in common events.
- The next and previous links may not correctly reflect the number of subject filtered in some instances. Going back to subject selection and re-selecting a patient solves the problem.
- Opening a form might fail in some instances when using IE 8.
- Empty CRF's PDF wont display all code list items for radio and checkboxes with many code list items in vertical layout.
- The code list tables for dropdowns are not displayed completely in the CRF PDF if the code list has many items.

Release Notes

Viedoc 4.9

Signoff



Signature Confirmation by
the product owner.

Uppsala, Sweden 2014-10-03

Henrik Blombergsson

PCG Solutions